DipSensor[®]

**RAPID TEST FOR THE DETECTION OF
BETA-LACTAMS (INCLUDING CEFALEXIN) AND
TETRACYCLINES IN COW'S MILK**

KIT900 | 25 ASSAYS

ENGLISH

FRANÇAIS

OPERATING PROCEDURE	3
HOW TO RUN A TEST	3
GENERAL REMARKS	5
RESULTS	6
VISUAL INTERPRETATION	6
TABLE A: PRELIMINARY ESTIMATE OF DETECTION LIMITS *	7

DISCLAIMER

The information contained in this document is provided in good faith, and every reasonable effort is made to ensure that it is correct and up to date.

However, given the dynamic nature of our business, and changing rules and regulations, there may be delays, omissions or inaccuracies in information contained in this document. Accordingly, Unisensor does not warrant the accuracy and completeness of the information.

The information contained in this document may be changed or updated at any time without notice. The information provided was last updated as of February 22, 2022.

Any person relying on any of the information contained in this document or making any use of the information contained herein, shall do so at its own risk.

To the best of our knowledge and belief, the products mentioned herein do not infringe the intellectual property rights of any third party. The products may be covered by pending or issued patents, registered or unregistered trademarks, or similar intellectual property rights.

OPERATING PROCEDURE

KIT900 is a rapid qualitative test that allows you to simultaneously detect and discriminate the presence of β -lactam and tetracycline molecules in a sample of raw cow's milk.

HOW TO RUN A TEST

1. Preparation


- Collect raw cow's milk in a container.
- Choose a clean and dry place to perform the test.
- Wash and dry your hands.

2. Immerse the Dipper tube in the milk, at least up to the half of the tube's height


3. Take the Dipper tube out of the milk container, and make sure there is milk in the bottom of the tube

4. Place the test strip in the Dipper tube in contact with the milk

OPERATING PROCEDURE

5. Incubate at room temperature for 10 minutes, to ensure sufficient migration of the reagents


6. Take the strip out of the Dipper tube and visually read results within 5 minutes
7. Wash the Dipper tube with plenty of water and leave to dry in a clean place
Do not reuse the tube in case of a positive result

GENERAL REMARKS

- The optimal environmental temperature range to run the test is 15 – 30 °C. If the environmental temperature is below 15°C, please add **5 additional minutes** to the test duration before conducting the interpretation. If under any circumstance the visual interpretation is doubtful, run the test for **15 minutes instead of 10 minutes**.
- Store the kit in the fridge at a temperature between 2°C and 8°C. Avoid repeated temperature variations.
- Before using the reagents, remove the kit from the refrigerator and wait for the reagents to reach room temperature.
- The tube must be clean before use.
- Avoid exposure of reagents to moisture and light.
- Do not use the kit after expiration date.
- Please avoid touching the membrane of the test strip. Always handle it from the top.
- Close the test strips container immediately after use.
- The cow's milk sample must be liquid and homogeneous. There can be no clots or sedimentation phases. The ideal temperature for the cow's milk sample is between 4°C and 37°C.
- Check that the strip reaches / touches the bottom of the dipper tube.
- Incubation time is 10 minutes. This timing must be respected to ensure a correct result.
- During the drying of the strip, the color intensities of the lines become sharper.
- When a positive result is recorded, that result should be confirmed by performing two more tests.

RESULTS

VISUAL INTERPRETATION

RESULTS

TABLE A: PRELIMINARY ESTIMATE OF DETECTION LIMITS *

B-LACTAMS (PPB)		TETRACYCLINES (PPB)	
PENICILLINES		Tetracycline	80 - 100
Ampicillin	2 - 4	Oxytetracycline	80 - 100
Penicillin G	2 - 3	Chlortetracycline	80 - 100
Amoxicillin	3 - 4	Doxycycline	80 - 100
Oxacillin	12 - 18		
Cloxacillin	6 - 8		
Dicloxacillin	9 - 10		
Nafcillin	80 - 110		
Piperacillin	1 - 4		
Ticarcillin	15 - 30		
Aspoxicillin	1 - 4		
Phenoxymethyl-penicillin (Pen V)	2 - 4		
Hetacillin	See Amoxi		
Penethamate	See PenG		
CEFALOSPORINES			
Ceftiofur	20 - 30		
Desfuroyl Ceftiofur	60 - 80		
Cefquinome	20 - 30		
Cefazolin	18 - 22		
Cefapirin	6 - 8		
Desacetyl Cephapirin	50 - 60		
Cefacetrile	30 - 40		
Cefoperazone	3 - 4		
Cefalexin	50 - 100		
Cefalonium	3 - 5		
Cefuroxime	140 - 200		

* These detection limit ranges should not be taken as definitive. They should be considered as indicative values for and they will be precisely defined and confirmed by subsequent validation studies.

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© Unisensor • All rights reserved